



LABORATORIES, INC.

286 ELDRIDGE ROAD, FAIRFIELD, NEW JERSEY 07004 U.S.A.

August 24, 2020

To Whom It May Concern:

Parker has performed functional and microbial testing to support heating the following products for 31 days at 43° C (or lower) in a Thermasonic® Gel Warmer:

Aquasonic® 100 Ultrasound Transmission Gel (01-08) 8 oz. bottles

Aquasonic® Clear Ultrasound Transmission Gel (03-08) 8 oz bottles

Scan® Ultrasound Transmission Gel (11-08) 8 oz bottles

Polysonic® Ultrasound Lotion (20-08) (21-08) 8.5 oz bottles

Aquasonic® 100 Ultrasound Transmission Gel Packet (01-20) 20 gram packettes

Sterile Aquasonic® 100 Ultrasound Transmission Gel Packet (01-01) 20 gram packettes

Aquasonic® Clear Ultrasound Transmission Gel Packet (03-20) 20 gram packettes

After 31 days, the gel product should be discarded.

Please let me know if I can be of any further assistance

Best regards,

A handwritten signature in purple ink, appearing to read "Adrina Lalabekova".

Adrina Lalabekova
Quality Assurance Specialist
Parker Laboratories, Inc.
alalabekova@parkerlabs.com